

# ÖSSUR FORMFIT® PRO KNEE FLITE

## INDICATIONS FOR USE

For knee conditions that may benefit from compression, such as:

- Mild sprains and strains
- Knee pain
- Feeling of instability

## FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit
- CoolVent™ knit in popliteal area for added breathability
- Patented FoldLock™ system for easy and secure application
- Lightweight, removeable spiral steel stays



## REQUIRED MEASUREMENTS



Measure calf circumference  
15 cm (6") below mid-patella.  
Leg in extension.



Measure thigh circumference  
15 cm (6") above mid-patella.  
Leg in extension.

Size	Calf circumference 15 cm (6") below mid patella		Thigh circumference 15 cm (6") above mid patella	
	cm	inch	cm	inch
XS	26 - 31	10 1/4 - 12 1/4	36 - 41	14 1/4 - 16 1/4
S	31 - 36	12 1/4 - 14 1/4	41 - 46	16 1/4 - 18
M	36 - 41	14 1/4 - 16 1/4	46 - 51	18 - 20
L	41 - 46	16 1/4 - 18	51 - 56	20 - 22
XL	46 - 50	18 - 19 3/4	56 - 60	22 - 23 1/2

## ORDERING INFORMATION

### BLUE

Size	Part#
XS	FP141001
S	FP141002
M	FP141003
L	FP141004
XL	FP141005



### BLACK

Size	Part#
XS	FP140001
S	FP140002
M	FP140003
L	FP140004
XL	FP140005